


Lifecare

SM_004 L

Capacity: -710mmHg \pm 10 at 32-35 LPM
Noise Level: <50 dB A \pm 3 Almost whispers
Pump type: Oil Free Double Piston Pump
Housing: Complete FRP Cabinet
Jars : Polycarbonate 2x2.5 liters capacity with mechanical over flow system & ABS Lid
Filter: Reusable bacterial filter
Tubing: Non - Collapsible Suction Tubing
Power: 220/230V AC, 50Hz, 1Ph
Vacuum Gauge: 2.0 inch, 0-760 mmHg

Rs. 16800/-




Instavac

SM_005 I

Capacity: -710mmHg \pm 10 at 30-35 LPM
Noise Level: <50 dB A \pm 3 Almost whispers
Pump type: Noiseless Rotary Vane Pump
Housing: Powder Coated MS Cabinet *
Jars : Polycarbonate 2x2.5 liters capacity with mechanical over flow system & ABS Lid
Filter: Reusable bacterial filter
Tubing: Non - Collapsible Suction Tubing
Power: 220/230V AC, 50Hz, 1Ph
Vacuum Gauge: 3.0 inch, 0-760 mmHg

Rs. 14700/-



Vacumaxx

SM_006 V

Capacity: -720mmHg \pm 10 at 32-35 LPM
Noise Level: <50 dB A \pm 3 Almost whispers
Pump type: Oil Free Double Piston Pump
Housing: Powder Coated Anti Corrosive MS Trolley
Jars : Polycarbonate 2x2.0 liters capacity with mechanical over flow system & ABS Lid
Filter: Reusable bacterial filter
Tubing: Non - Collapsible Suction Tubing
Power: 220/230V AC, 50Hz, 1Ph
Vacuum Gauge: 2.0 inch, 0-760 mmHg

Rs. 21000/-



* Carrying Warranty at Rajkot Office only on site warranty not provided free of cost.

* Due to constant upgradation design, price and features are subject to change any time without prior notice.



Western Surgical

WESTERN SURGICAL

www.westernsurgical.in

Email : westernsurgical@gmail.com

Ph. : +91-281-2232871 / 2243159

2

© Copyright Reserved by : Western Surgical - Rajkot.