

Gynaecology colposcope P 34 A

Green Filter Helps to
view Blood Vessels/
Tissues/ skin colours
without discolouration

Cold Light Source
Halogen Fibre - Optic Light
Source. No heat/warmth to
inconvenience patient or
examiner. Operates on
230V AC. Available with
Single or Double outlet.

**Exclusive Tilting
arrangement**
Viewing angle can be
changed to suit seated,
standing or
stooping viewers/
Surgeons

Note:
Colposcope
is supplied
duly fitted
with the stan-
dard 300 mm
Objective Lens.
Other two Lenses
(200mm & 400mm) will
be supplied as optional at
extra cost. Colposcope is
supplied with single
outlet light source. Twin
light source will be
optional at extra cost.

Model CM-51020