

TEAR OPEN HERE

SURGICARE®

STERILE

**DISPOSABLE
SURGICAL
RUBBER GLOVES**
PRE-POWDERED
MICRO ROUGH TEXTURE



CE 2400
0020

ONE PAIR (2N)
③ SINGLE USE



CAUTION:

After donning remove powder by wiping
gloves thoroughly with a sterile wet
sponge or sterile wet cloth or other
appropriate method.
This product contains natural rubber
latex which may cause a latex reaction,
including anaphylaxis, in some individuals.
Do not use if package is puffed
damaged or wet.
*Do not reuse. Single use only.
Check for latex and allergen content.

kanam

KANAM LATEX INDUSTRIES PRIVATE LIMITED
ISO 9001:2008, EN ISO 13485:2012 and ISO 14001:2004 Certified Company