

TEAR HERE TO OPEN

# SURGICARE®



TYPE 1  
GM/L 8105290



0334

**DISPOSABLE  
SURGICAL  
RUBBER GLOVES**  
PPE-POWDERED  
MICRO ROUGH TEXTURE

ONE PAIR (2 Nos)

SINGLE USE

**STERILE**

**SIZE : S**

Width : 102 ± 6mm

**CAUTION**

After donning remove powder by  
wiping gloves thoroughly with  
sterile alcohol prep pad or saline wet  
cloth or other effective method.

This product contains natural rubber  
latex which may cause allergic  
reactions, including anaphylactic  
reactions to some individuals.

Do not use if packaging is  
perforated, damaged, or aged.

Do not reuse. Reuse can  
cause cross-infection and  
compromise safety.



**KANAM LATEX INDUSTRIES PRIVATE LIMITED**

ISO 9001:2008, ISO 13485:2003 and ISO 14001:2014 Certified Company

M