

# Endo / Laparoscopy Trolley



## Special Features

1. Size - 4.75 Ft. H x L + W - 2.25 Feet
2. Attractive Color Combination
3. Self Size - 18" x 24" tray type
4. Self Height Adjustable No. 2,3,4,5
5. One Drawer with Lock
6. Stand for Cautry foot switch
7. 4" Castor with Lock systems
8. Two S.S. Handle for Movement
9. Co2 Cylinder stand Back Side with 2 clamp
10. 6 + 1 Plug 15 x 5 Amp Switch & one  
Extra for Cautry with 10 Ft. Cable
11. Space for UPS / Servo Stabilizer behind



## DUAL CABLE CONNECTION

18,000/-

\* Carrying Warranty at Rajkot Office only on site warranty not provided free of cost.

\* Due to constant upgradation design, price and features are subject to change any time without prior notice.



**WESTERN SURGICAL**  
EQUIPMENTS PRIVATE LIMITED

Western Surgical

\* [www.westernsurgical.in](http://www.westernsurgical.in) \*

Email : [westernsurgical@gmail.com](mailto:westernsurgical@gmail.com)

Ph. +91-281-2232871 / 2243159

Mo. +91 94280 03750/52/54

© Copyright Reserved by : Western Surgical - Rajkot.