

Rao's Hand Suction Unit 100 ml Capacity

1102



Rao's Hand Suction Unit 40 ml Capacity

1101



The Hand Suction Units are suitable for mucous extraction or to supplement general resuscitation techniques where comprehensive suction facilities are not available. The containers are unbreakable and transparent, permitting immediate visual examination of the aspirate.

CE



Goldman Halothane Vaporizer

4002

Without stand & adaptors - for use only with Anaesthesia Apparatus.

4001

With stand, adaptors & oxygen tubing with plug-in mount for general use.

The Goldman Halothane Vaporizer is designed to limit the maximum concentration of Halothane Vapour to 3% by volume, irrespective of the total gas flow through the unit. Adjustment is made by means of a ratchet type control valve. The Glass Jar is screwed onto the body and can be quickly detached for easy filling. The Vaporizer is also suitable for use with Trichloroethylene in low concentrations.

* Carrying Warranty at Rajkot Office only on site warranty not provided free of cost.

* Due to constant upgradation design, price and features are subject to change any time without prior notice.



Western Surgical

WESTERN SURGICAL
EQUIPMENTS PRIVATE LIMITED

* www.westernsurgical.in *

Email : westernsurgical@gmail.com

Ph. +91-281-2232871 / 2243159

Mo. +91 94280 03750/52/54

© Copyright Reserved by : Western Surgical - Rajkot.