


# AVSU ZONAL VALVE BOX

Zone Valves and Zone Valve Boxes provide versatility and value through medical gas or vacuum pipeline flow prevention.



## Features

- Offered in single (houses 1 valve) or multiple (houses 2 to 6 valves) types.
- Features a pull-out, removable opaque window.
- Adjustable to different wall thicknesses.
- Include color-coded gas identification labels.
- Blowout proof stem.
- Available in locking or non-locking.
- Swing-out design to prevent damage during brazing operation.
- 100% hydrostatically tested.
- 3-piece ball type design for ease of maintenance.

## MANIFOLD SYSTEM (OXYGEN/N2O)

The manifolds manufactured by us are extendable & Non-extendable type with choice of size of cylinders bank for bulk 'D' type jumbo cylinders.

Manifold is suitable to withstand a working pressure of 140Kg/cm<sup>3</sup>.

Manifolds have high pressure copper annealed tail pipes with one end having Brass adaptor suitable for cylinders. Manifolds have Brass Non Return Valves (NRV) for individual cylinders (Cylinders to be provided by the client).



\* Carrying Warranty at Rajkot Office only on site warranty not provided free of cost.

\* Due to constant upgradation design, price and features are subject to change any time without prior notice.



Western Surgical

**WESTERN SURGICAL**  
EQUIPMENTS PRIVATE LIMITED

45

\* [www.westernsurgical.in](http://www.westernsurgical.in) \*  
Email : [westernsurgical@gmail.com](mailto:westernsurgical@gmail.com)  
Ph. +91-281-2232871 / 2243159  
Mo. +91 94280 03750/52/54

© Copyright Reserved by : Western Surgical - Rajkot.