

Heine Otoscopes Full Size Instruments P 415 A



HEINE BETA 200[®] Otoscope
Distal Fiber optics (F.O.) illumination.
Metal body with 3 x magnifying window.
Page 6



HEINE K 180[®] Otoscope
As BETA 200 with distal F. O. illumina-
tion. High impact polycarbonate body.
The economical choice.
Page 6



**HEINE BETA[®] 100 Diagnostic
Otoscope**
A new, all-metal otoscope with direct ha-
logen illumination. Unusually versatile in
combination with nasal- and spreadable
specula for easy instrumentation. Matt
black interior eliminates reflexes.



**HEINE K 100[®] Diagnostic
Otoscope**
A new, competitively-priced otoscope
with direct halogen illumination. High-
impact black polycarbonate body with a
connector for nasal- and spreadable
specula for easy instrumentation.



HEINE Operating Otoscope
Universal instrument – ideal for
instrumentation under magnification.
The loupe and speculum can be
swivelled to one side. Not suitable for
pneumatic testing. Page 7