



Powervac

SM_007 P

Capacity: -730mmHg $\pm$ 10 at 60 LPM
Noise Level: <50 dB A $\pm$ 3 Almost whispers
Pump type: Oil Free Double Piston Pump
Housing: Powder Coated MS Cabinet
Jars : Polycarbonate 2x2.5 liters capacity with mechanical over flow system & ABS Lid
Filter: Reusable bacterial filter
Tubing: Non - Collapsible Suction Tubing
Power: 220/230V AC, 50Hz, 1Ph
Vacuum Gauge: 2.0 inch, 0-760

M.R.P. 60000/-

Smart Lipovac

SM_008 SL

Capacity: -740mmHg $\pm$ 10 at 120 LPM
Noise Level: <50 dB A $\pm$ 3 Almost whispers
Pump type: Oil Free Double Piston Pump
Housing: Powder Coated MS Cabinet
Jars : Polycarbonate 2x2.5 liters capacity with mechanical over flow system & ABS Lid
Filter: Reusable bacterial filter
Tubing: Non - Collapsible Suction Tubing
Power: 220/230V AC, 50Hz, 1Ph
Vacuum Gauge: 4.0 inch, 0-760 mmHg

M.R.P. 75000/-



AirForce

Housing : Powder Coated Anti Corrosive MS
Jars : Polycarbonate 2 X 4.0 liters capacity
Overflow Protection : Mechanical Type
Tubing : Nom - Collapsible Suction Tubing
Vacuum Gauge : 2.5 inch, 0-760 mmHg
Filter : Bacterial Filter Autoclavable / Reusable

Pump Type : Oil Free Twin Headed Piston Pump
Capacity : 740mmHg $\pm$ 10 at 120~125 LPM
Noise Level : <50 dB A $\pm$ 3 Almost whispers
Motor : 220/230V AC, 50 Hz, 820W
Valves : Stainless Steel

M.R.P. 115000/-



* Carrying Warranty at Rajkot Office only on site warranty not provided free of cost.

* Due to constant upgradation design, price and features are subject to change any time without prior notice.



Western Surgical

WESTERN SURGICAL
EQUIPMENTS PRIVATE LIMITED

* www.westernsurgical.in *

Email : westernsurgical@gmail.com

Ph. +91-281-2232871 / 2243159

Mo. +91 94280 03750/52/54

© Copyright Reserved by : Western Surgical - Rajkot.