



Imported China

SM_0028 IMC

- + Adopting oil free piston pump
- + Clean without pollution of oil and smoke
- + Convent no need add oil to maintain daily
- + Safe no plus pressure during using
- + No flow backwards pressure when the machine stops.
So the liquid won't flow backwards
- + Complete plastic panels, hand-switch and foot-switch
for operating easily. Select negative pressure at will
with no level for pressure adjustment, and overflow
protection mechanism, safe and reliable
- + with trolley type folding handle

Principle Technical Data

max negative pressure $\geq 0.09\text{MPa}$

noise $\leq 65\text{dB (A)}$

Power $\leq 180\text{VA}$

pumping rate $\geq 20\text{L / min}$

reservoir capacity 2500mL / pc. 2 pieces

Note

not suitable for continuous use in the induced abortions
and for the inflammable or explode surroundings. Max
continuous working time in 30 minutes, continuity rate 50%

electrical appliance safety requirement class I, type B

M.R.P. 30000/-

* Carrying Warranty at Rajkot Office only on site warranty not provided free of cost.

* Due to constant upgradation design, price and features are subject to change any time without prior notice.



Western Surgical

WESTERN SURGICAL
EQUIPMENTS PRIVATE LIMITED

* www.westernsurgical.in *

Email : westernsurgical@gmail.com

Ph. +91-281-2232871 / 2243159

Mo. +91 94280 03750/52/54

© Copyright Reserved by : Western Surgical - Rajkot.