


ONE STOP SHOWROOM FOR HOSPITAL & HEALTHCARE PRODUCTS.



SUCTION MACHINE

CATALOGUE - 2023



Western[®] Surgical

WESTERN SURGICAL[®]

AN ISO 9001 : 2008 & ISO 13485 : 2015 AC 2012, CE, FDA & GMP CERTIFIED COMPANY

www.westernsurgical.in Email : westernsurgical@gmail.com

Our Valuable Clients

Sr. No.	Dr. Name & Hospital Name	City Name
1	Dr. Viral Vasavada (general Surgeon) M. S. Samarpan Hospital	RAJKOT
2	Dr. Sanjay Pandya M.D., D.N.B., Nephro	RAJKOT
3	Dr. Raju Sagar M.S. Laparoscopy Surgeon	RAJKOT
4	Dr. I.P. Swami M.D. Physician	RAJKOT
5	Dr. Pradip Shah Pediatric	RAJKOT
6	Dr. Rajesh Gandhi Orthopedic	RAJKOT
7	Dr. Rakesh Chokshi Laproscopy Surgeon	RAJKOT
8	Dr. Sanjay Patoliya M.S. Bariatric Surgeon	RAJKOT
9	Dr. Shrenik Shah M.S. ENT	RAJKOT
10	Dr. Himanshu Thakkar M.S. ENT	RAJKOT
11	Dr. Samir Thakkar Pediatric	RAJKOT
12	Dr. Vinay Majmudar M.D. Physician	RAJKOT
13	Dr. Hitesh Dhru M.S. MCH. Plastic	RAJKOT
14	Dr. Girish Amlani M.S. MCH. Plastic	RAJKOT
15	Dr. Bhaumik Bhayani M.S. MCH. Plastic	RAJKOT
16	Dr. Ashish Jasani M.S. MCH. Urologist	RAJKOT
17	Dr. Bhavesh Sachde Orthopedic Join Replacement Surgeon	RAJKOT
18	Guru Govindsinh Hospital General Hospital	JAMNAGAR
19	Civil Hospital - RAJKOT, AMRELI, SURAT, AHMEDABAD, MORBI, PORBANDAR, ETC.	
20	Aashirvad Multispicility Hospital	RAJKOT
21	Bhavsinhji General Hospital	PORBANDAR
22	Sir-t Hospital General Hospital	BHAVNAGAR
23	Swami Nirdoshanand Hospital - Timbi Multispeciality Hospital	TIMBI
24	Samarpan General Hospital	JAMNAGAR
25	C. U. Shah Medical College	SURENDRANAGAR
26	Christ Hospital General Hospital	RAJKOT
27	Sadbhavna Hospital	MAHUVA
28	City Dental Hospital	RAJKOT
29	Giriraj Multispeciality Hospital	RAJKOT
30	H. J. Doshi Sarvajanic Hospital General Hospital	RAJKOT

For More Valuable Customer on our web site www.westernsurgical.in

Our Valuable Clients

Sr. No.	Dr. Name & Hospital Name	City Name
31	Hanumant Hospital Multispeciality Hospital	MAHUVA
32	Shri Ranchoddas Babu Aasharm Eyes Hospital	RAJKOT
33	Wockhardt Hospital General Hospital	RAJKOT
34	B. A. Danger Homeopethic College Medical College	RAJKOT
35	Dr. Nilang Vasavda M.D. Medicine	RAJKOT
36	Dr. Hemang Vasavda & Dr. Hurd Vasavda	RAJKOT
37	Dr. Punit Trivedi M.S. MOU	RAJKOT
38	Dr. Avinash Maru M.S. Ortho.	RAJKOT
39	Dr. Gaurang Dave M.S. Surgeon	RAJKOT
40	Krishna Multispeciality Hospital	MORBI
41	Ayush Multispeciality Hospital	MORBI
42	Dr. Janak Thakkar M.D. Skin	RAJKOT
43	The Gujarat Cancer & Reserch Institute Research Institute	AHMEDABAD
44	U. N. Mehta Institue Of Cardiology And Research Centre Research Institute	AHMEDABAD
45	Sterling Hospital Multispeciality Hospital	RAJKOT, GANDHINAGAR
46	HCG Hospital Multispeciality Hospital	RAJKOT
47	Olympus Hospital Multispeciality Hospital	RAJKOT
48	Synergy Superspeciality Hospital Multispeciality Hospital	RAJKOT
49	B. T. Savani Kidney Hospital Kidney Hospital	RAJKOT
50	Kakrapar Atomic Power Station Hospital General Hospital	KAKRAPAR, VYARA
51	Shree G. T. Sheth Orthopaedic Hospital Joint Replac. Surgen	RAJKOT
52	Sadbhavna Hospital Multispeciality Hospital	RAJKOT
53	Dasha Shrimali Jain Hospital General Hospital	RAJKOT
54	Dr. Vinod Kaila M.D. Physician	MORBI
55	Dr. Dhaval H. Mankad M.S. Surgeon	JAMNAGAR
56	Shree Jalaram Raghukul Sarvajanik Multispeciality Hospital	RAJKOT
57	Orthoplus Speciality Hospital Orthopedic Joint Replac. Surgen	RAJKOT
58	Dr. Nita Thakkar Gaynec Aagman Maternity Home	RAJKOT
59	Dr. Harshal Mistri	SURAT

For More Valuable Customer on our web site www.westernsurgical.in

INDEX

Minic, Nano - AC/DC, Medevac.....	1
Lifecare, Instavac, Vacuumaxx.....	2
Powervac, Smart Lipovac.....	3
Supreme AirForce.....	3
Foot Suction Reg. S.S. Top, Foot Suction Imported.....	4
Vacuum Extractor, Silicone Cup.....	4
Principle Technical Data.....	5
Accessories for Suction Machine.....	6
Accessories for Suction Machine.....	7
Suction Machine S.S. Top.....	8
Suction Machine Fully Steel Body.....	8
Suction Foot Switch.....	9

* Carrying Warranty at Rajkot Office only on site warranty not provided free of cost.

* Due to constant upgradation design, price and features are subject to change any time without prior notice.



WESTERN SURGICAL

Western Surgical

www.westernsurgical.in

Email : westernsurgical@gmail.com

Ph. : +91-281-2232871 / 2243159

© Copyright Reserved by : Western Surgical - Rajkot.



Minic SM_001 M

Capacity: -550 mmHg \pm 10 at 12-15 LPM
Noise Level: <50 dB A \pm 3 Almost whispers
Pump type: Oil Free Piston Pump
Housing: Shock Proot FRP Cabinet
Jars : Polycarbonate 2x2.5 liters capacity with mechanical over flow system & ABS Lid
Filter: Reusable bacterial filter
Tubing: Non - Collapsible Suction Tubing
Power: 220/230V AC, 50Hz, 1Ph
Vacuum Gauge: 2.0 inch, 0-760 mmHg

M.R.P. 12500/-

Nano - AC/DC

SM_002 N AC/DC

Capacity: -580 mmHg \pm 10 at 12-15 LPM
Backup : 120 mins
Pump type: Oil Free Piston Pump
Housing: Complete ABS (Plastic) Cabinet
Jars : Polycarbonate 2x2.0 liters capacity with mechanical over flow system & ABS Lid
Filter: Reusable bacterial filter
Tubing: Non - Collapsible Suction Tubing
Power: 220/230V AC, 50Hz, 1Ph & 12V DC
Vacuum Gauge: 2.0 inch, 0-760 mmHg

M.R.P. 21000/-



Medevac

SM_003 M

Capacity: -710mmHg \pm 10 at 28-32 LPM
Noise Level: <50 dB A \pm 3 Almost whispers
Pump type: Noiseless Rotary Vane Pump
Housing: Powder Coated MS Cabinet
Jars : Glass 2x2.0 liters capacity with mechanical over flow system & ABS Lid
Filter: Reusable bacterial filter
Tubing: Non - Collapsible Suction Tubing
Power: 220/230V AC, 50Hz, 1Ph
Vacuum Gauge: 2.5 inch, 0-760 mmHg

M.R.P. 23500/-



* Carrying Warranty at Rajkot Office only on site warranty not provided free of cost.
 * Due to constant upgradation design, price and features are subject to change any time without prior notice.



WESTERN SURGICAL

www.westernsurgical.in
 Email : westernsurgical@gmail.com
 Ph. : +91-281-2232871 / 2243159



© Copyright Reserved by : Western Surgical - Rajkot.

Instavac

SM_005 I

Capacity: -710mmHg \pm 10 at 30-35 LPM

Noise Level: <50 dB A \pm 3 Almost whispers

Pump type: Noiseless Rotary Vane Pump

Housing: Powder Coated MS Cabinet *

Jars : Polycarbonate 2x2.5 liters capacity with mechanical over flow system & ABS Lid

Filter: Reusable bacterial filter

Tubing: Non - Collapsible Suction Tubing

Power: 220/230V AC, 50Hz, 1Ph

Vacuum Gauge: 3.0 inch, 0-760 mmHg

M.R.P. 50000/-



Vacumaxx

SM_006 V

Capacity: -720mmHg \pm 10 at 32-35 LPM

Noise Level: <50 dB A \pm 3 Almost whispers

Pump type: Oil Free Double Piston Pump

Housing: Powder Coated Anti Corrosive MS Trolley

Jars : Polycarbonate 2x2.0 liters capacity with mechanical over flow system & ABS Lid

Filter: Reusable bacterial filter

Tubing: Non - Collapsible Suction Tubing

Power: 220/230V AC, 50Hz, 1Ph

Vacuum Gauge: 2.0 inch, 0-760 mmHg

M.R.P. 55000/-

* Carrying Warranty at Rajkot Office only on site warranty not provided free of cost.

* Due to constant upgradation design, price and features are subject to change any time without prior notice.



WESTERN SURGICAL

Western Surgical


www.westernsurgical.in

Email : westernsurgical@gmail.com

Ph. : +91-281-2232871 / 2243159

2

© Copyright Reserved by : Western Surgical - Rajkot.



Powervac

SM_007 P

Capacity: -730mmHg \pm 10 at 60 LPM
Noise Level: <50 dB A \pm 3 Almost whispers
Pump type: Oil Free Double Piston Pump
Housing: Powder Coated MS Cabinet
Jars : Polycarbonate 2x2.5 liters capacity with mechanical over flow system & ABS Lid
Filter: Reusable bacterial filter
Tubing: Non - Collapsible Suction Tubing
Power: 220/230V AC, 50Hz, 1Ph
Vacuum Gauge: 2.0 inch, 0-760

M.R.P. 60000/-

Smart Lipovac

SM_008 SL

Capacity: -740mmHg \pm 10 at 120 LPM
Noise Level: <50 dB A \pm 3 Almost whispers
Pump type: Oil Free Double Piston Pump
Housing: Powder Coated MS Cabinet
Jars : Polycarbonate 2x2.5 liters capacity with mechanical over flow system & ABS Lid
Filter: Reusable bacterial filter
Tubing: Non - Collapsible Suction Tubing
Power: 220/230V AC, 50Hz, 1Ph
Vacuum Gauge: 4.0 inch, 0-760 mmHg

M.R.P. 75000/-



AirForce

Housing : Powder Coated Anti Corrosive MS
Jars : Polycarbonate 2 X 4.0 liters capacity
Overflow Protection : Mechanical Type
Tubing : Nom - Collapsible Suction Tubing
Vacuum Gauge : 2.5 inch, 0-760 mmHg
Filter : Bacterial Filter Autoclavable / Reusable

Pump Type : Oil Free Twin Headed Piston Pump
Capacity : 740mmHg \pm 10 at 120~125 LPM
Noise Level : <50 dB A \pm 3 Almost whispers
Motor : 220/230V AC, 50 Hz, 820W
Valves : Stainless Steel

M.R.P. 115000/-

* Carrying Warranty at Rajkot Office only on site warranty not provided free of cost.
* Due to constant upgradation design, price and features are subject to change any time without prior notice.



WESTERN SURGICAL

Western Surgical

www.westernsurgical.in

Email : westernsurgical@gmail.com

Ph. : +91-281-2232871 / 2243159

3

© Copyright Reserved by : Western Surgical - Rajkot.



Foot Suction Regular S.S. Top

SM_009 FSRT

M.R.P. 6000/-

Portable Twin glass of 1 Litr superior piston pump for creating vacuum instantly, base made from M.S. Sheet full powder coated, all aluminum parts are fully powder coated cacuum guage 2½" guage 2½" graduated in mmHg. (Also available in stainless steel base.)



Foot Suction Imported

SM_0011 FSI

Specification

Base : Power Coated Heavy MS Base
 Capacity : -600mg Hg±10 at 200ml Per Stroke
 Pump Type : Heavy Piston Pump
 Jar : Polycarbonate 1500ml capacity
 : with over flow system & ABS Lid
 Filter : Bacterial Filter Autoclavable/Reusable
 Vacuum Gauge : 2.0 inch, 0-760mm Hg
 Operated : Manual

M.R.P. 9500/-



Vacuum Extractor

SM_0012 VE

M.R.P. 15000/-

- Supreme Manual Vacuum Extractor
- Traction handle & Hand Suction pump
- Three Silicone cups (40,50,60mm)
- Glass jar fitted with vacuum gauge and control valve
- Optional:Stainless Steel instead of Silicone cups can be provided.
- Optional:Rubber tubes instead of silicone tubes can be provided.
- Stand for glass jar / • Silicone pressure tube



Silicone Cup SM_0013 SC

- 100% silicone medical product
- This vacuum suction cup could assist and shorten second stage of labor in occipital delivery
- Minimize the risk of Traumas, Cephalohematoma and Forceps delivery
- The soft cup can protect and decrease possible danger of damaging the fetus' head
- The handle and valve assemblies are manufactured of chrome plated brass and can be replaced

M.R.P. 3500/-

* Carrying Warranty at Rajkot Office only on site warranty not provided free of cost.

* Due to constant upgradation design, price and features are subject to change any time without prior notice.



WESTERN SURGICAL

Western Surgical

www.westernsurgical.in

Email : westernsurgical@gmail.com

Ph. : +91-281-2232871 / 2243159

4

© Copyright Reserved by : Western Surgical - Rajkot.



Imported China

SM_0028 IMC

- + Adopting oil free piston pump
- + Clean without pollution of oil and smoke
- + Convenient no need add oil to maintain daily
- + Safe no plus pressure during using
- + No flow backwards pressure when the machine stops.
So the liquid won't flow backwards
- + Complete plastic panels, hand-switch and foot-switch for operating easily. Select negative pressure at will with no level for pressure adjustment, and overflow protection mechanism, safe and reliable
- + with trolley type folding handle

Principle Technical Data

max negative pressure $\geq 0.09\text{MPa}$
noise $\leq 65\text{dB (A)}$
Power $\leq 180\text{VA}$
pumping rate $\geq 20\text{L / min}$
reservoir capacity 2500mL / pc. 2 pieces

Note

not suitable for continuous use in the induced abortions and for the inflammable or explode surroundings. Max continuous working time in 30 minutes, continuity rate 50%
electrical appliance safety requirement class I, type B

M.R.P. 30000/-

* Carrying Warranty at Rajkot Office only on site warranty not provided free of cost.
 * Due to constant upgradation design, price and features are subject to change any time without prior notice.



WESTERN SURGICAL

Western Surgical

www.westernsurgical.in

Email : westernsurgical@gmail.com

Ph. : +91-281-2232871 / 2243159

5

© Copyright Reserved by : Western Surgical - Rajkot.

Accessories for suction Machine

SM_0029 BF



Rs. 210/-

Bacterial Filter

SM_0030 PJ



Rs. 2100/-

2000ml PC Jar

SM_0031 HF



Rs. 315/-

Hydrophobic Filter

SM_0032 VP



M.R.P. 6000/-

Vacuum Pump



Rs. 630/-

SM_0033 FS

Foot Switch

SM_0034 TCP



Rs. 210/-

Twin Castor PU

* Carrying Warranty at Rajkot Office only on site warranty not provided free of cost.
* Due to constant upgradation design, price and features are subject to change any time without prior notice.



WESTERN SURGICAL

Western Surgical

www.westernsurgical.in
Email : westernsurgical@gmail.com
Ph. : +91-281-2232871 / 2243159

Accessories for suction Machine

SM_0035 PJ



SM_0036 CSP

Rs. 315/-



Cap For Suction Bottel PVC

SM_0037 CSG

Rs. 315/-



Cap For Suction Bottel Glass

SM_0038 CSG

Rs. 210/-



Cap For Suction Bottel Glass



Rs. 1050 to 2100/-



Polycarbonate Jars



SM_0041 GJ

Rs. 840/-

2500ml Glass Jar



SM_0042 PVC J

Rs. 1050/-

600ml PVC Jar



SM_0043 GJ

Rs. 840/-

1500ml Glass Jar

* Carrying Warranty at Rajkot Office only on site warranty not provided free of cost.
 * Due to constant upgradation design, price and features are subject to change any time without prior notice.



WESTERN SURGICAL

Western Surgical

www.westernsurgical.in

Email : westernsurgical@gmail.com

Ph. : +91-281-2232871 / 2243159

7

© Copyright Reserved by : Western Surgical - Rajkot.

SM_0039 SED

Features

- Superior Quality Suction Pump for Instant Vacuum
- Unique removable filter
- Heavy duty anti static castors
- Non collapsible 10mm tube
- Versatile fo M.T.P. & othe major surgical applications
- Oil rotary noseless pump

Specification

Vacuum	: -700+10mm Hg control by metal knob
Sound	: 50dB A+B almost noiseless
Suction Capacity	: 35Ltr/min
Jars	: 2x2Ltrs Glass jars with over flow safety Device
Body	: M.S. Body S.S. Top
Gauge	: 3" graduated in mm Hg
Voltage	: 220V 50 C/s
Motor	: ½ HP

- * Also available Poly Carbonate Jar
- * Available in Crompton Motor & Local Motor

M.R.P. 30000/-



Suction Machine Fully Steel Body Oil Emergen Pump Pully Drive

SM_0040 SBED

Features

- Superior Quality Suction Pump for Instant Vacuum
- Heavy duty anti static castors
- Versatile fo M.T.P. & other major surgical applications
- Unique removable filter
- Non collapsible 10mm tube
- Oil rotary noseless pump

Specification

Vacuum	: -700+10mm Hg control by metal knob
Sound	: 50dB A+B almost noiseless
Suction Capacity	: 35Ltr/min
Jars	: 2x2Ltrs Glass jars with over flow safety Device
Body	: Fully Steel
Gauge	: 3" graduated in mm Hg
Voltage	: 220V 50 C/s
Motor	: ¼ HP

- * Also available Poly Carbonate Jar
- * Available in Crompton Motor & Local Motor

M.R.P. 45000/-



* Carrying Warranty at Rajkot Office only on site warranty not provided free of cost.
* Due to constant upgradation design, price and features are subject to change any time without prior notice.



WESTERN SURGICAL

Western Surgical

www.westernsurgical.in

Email : westernsurgical@gmail.com

Ph. : +91-281-2232871 / 2243159

8

© Copyright Reserved by : Western Surgical - Rajkot.

SUCTION FOOT SWITCH



Rs. 360/-

SUCTION FOOT SWITCH



Rs. 785/-

* Carrying Warranty at Rajkot Office only on site warranty not provided free of cost.
* Due to constant upgradation design, price and features are subject to change any time without prior notice.





WESTERN SURGICAL

Western Surgical

www.westernsurgical.in
Email : westernsurgical@gmail.com
Ph. : +91-281-2232871 / 2243159

OUR OTHER PRODUCTS RANGE



* Carrying Warranty at Rajkot Office only on site warranty not provided free of cost.

* Due to constant upgradation design, price and features are subject to change any time without prior notice.



WESTERN SURGICAL

Western Surgical

www.westernsurgical.in

Email : westernsurgical@gmail.com

Ph. : +91-281-2232871 / 2243159

© Copyright Reserved by : Western Surgical - Rajkot.

OUR OTHER PRODUCTS RANGE



* Carrying Warranty at Rajkot Office only on site warranty not provided free of cost.

* Due to constant upgradation design, price and features are subject to change any time without prior notice.

Dt. 02-12-2023

Printing Dt. 02-12-2023



WESTERN SURGICAL

Western Surgical

www.westernsurgical.in

Email : westernsurgical@gmail.com

Ph. : +91-281-2232871 / 2243159

© Copyright Reserved by : Western Surgical - Rajkot.



Western[®] Surgical

WESTERN SURGICAL[®]

We are Available on :        GeM
Government
e Marketplace

AN ISO 9001 : 2008 & ISO 13485 : 2015 AC 2012, CE, FDA & GMP CERTIFIED COMPANY

"Happy Home", 1st Floor, 23,
Karanpara, B/h. S.T. Bus Stand,
Kanak Road, Rajkot-360001.

Ph. : +91-281-2232871/2243159,

Mo. : +91 98250 69069 - +91 98250 45997

Email : westernsurgical@gmail.com

Email : info@westernsurgical.in

* www.westernsurgical.in *

© Copyright Reserved by : Western Surgical - Rajkot.



OUR MANUFACTURING & SHOWROOM UNIT

