

Suction Machine S.S. Top Oil Emergen Pump Pully Drive

SM_0039 SED

Features

- Superior Quality Suction Pump for Instant Vacuum
- Unique removable filter
- Heavy duty anti static castors
- Non collapsible 10mm tube
- Versatile fo M.T.P. & othe major surgical applications
- Oil rotary noseless pump

Specification

Vacuum : -700+10mm Hg control by metal knob
Sound : 50dB A+B almost noiseless
Suction Capacity : 35Ltr/min
Jars : 2x2Ltrs Glass jars with over flow safety Device
Body : M.S. Body S.S. Top
Gauge : 3" graduated in mm Hg
Voltage : 220V 50 C/s
Motor : 1/2 HP

- * Also available Poly Carbonate Jar
- * Available in Crompton Motor & Local Motor

M.R.P. 30000/-



Suction Machine Fully Steel Body Oil Emergen Pump Pully Drive

SM_0040 SBED

Features

- Superior Quality Suction Pump for Instant Vacuum
- Heavy duty anti static castors
- Versatile fo M.T.P. & other major surgical applications
- Unique removable filter
- Non collapsible 10mm tube
- Oil rotary noseless pump

Specification

Vacuum : -700+10mm Hg control by metal knob
Sound : 50dB A+B almost noiseless
Suction Capacity : 35Ltr/min
Jars : 2x2Ltrs Glass jars with over flow safety Device
Body : Fully Steel
Gauge : 3" graduated in mm Hg
Voltage : 220V 50 C/s
Motor : 1/4 HP

- * Also available Poly Carbonate Jar
- * Available in Crompton Motor & Local Motor

M.R.P. 45000/-



* Carrying Warranty at Rajkot Office only on site warranty not provided free of cost.

* Due to constant upgradation design, price and features are subject to change any time without prior notice.



Western Surgical

WESTERN SURGICAL
EQUIPMENTS PRIVATE LIMITED

* www.westernsurgical.in *

Email : westernsurgical@gmail.com

Ph. +91-281-2232871 / 2243159

Mo. +91 94280 03750/52/54

© Copyright Reserved by : Western Surgical - Rajkot.