

WS - 203 TELESCOPIC DELIVERY TABLE

Size :72"x 30" x 30" ht.



Tubular frame work powder coated or fully S.S. with 18 Swg. stainless steel Top. Provided with saline rod, thigh rest, stainless steel bowl. Telescopic foot end on smooth heavy duty sliding channels, Backrest adjustable on ratchet and head low with gas spring mechanism. Provided with powder coated or S.S. Detachable side rails and head low.

Attendance Bed with Locker Table



6000/-

WSAB2L : Attendance Bed & 2 Locker




7600/-

WSAB-4LWM : Attendance Bed - 4L Without Mattress



8500/-

WSAB4LM : Attendance Bed & Mattress



10,500/-

* Carrying Warranty at Rajkot Office only on site warranty not provided free of cost.

* Due to constant upgradation design, price and features are subject to change any time without prior notice.



WESTERN SURGICAL
EQUIPMENTS PRIVATE LIMITED

Western Surgical

36

* www.westernsurgical.in *

Email : westernsurgical@gmail.com

Ph. +91-281-2232871 / 2243159

Mo. +91 94280 03750/52/54

© Copyright Reserved by : Western Surgical - Rajkot.