

**WSP-111-D
Single Arm
Surgical / Anaesthesia
Pendant**



Single Arm Surgical Pendant Is New Addition In Field Medical Equipment In Modern Operating Theater : This Medical Pendant Has A Compact Structure By Logical Design. The Arm Console Head Can Be Rotated 340° The Pendant Is Manufactured In Accordance To The Ceiling Height Of The Operation Room And In Consideration With The False Ceiling To Facilitate Easy Movement Of The Pendant For The Medical Staff. Provided With 2 Shelves and Drawer System.

Some Necessary Medical Equip. Like Medical Gas Outlets And High/low Voltage Sockets. Internet Ports Are Integral Part In Medical Pendants. Also Provide IV Pole Attachments And Baskets In SS



**WSP-112-E
Wall / Ceiling Mounted
I.C.U. Pendant**



- 1) Provision For Gas Outlets
2 Oxygen 1 Vacuum, 1 Air.
- 2) 6 Electrical Connections.
- 3) 1 RJ 45
- 4) 4 Pole IV Attachment.
- 5) Monitor Arm With Monitor Holder.
- 6) 4 Shelves For Syringe/ Infusion Pump.
- 7) Wall / Ceiling Mounted

Extra Accessories Like Gas Outlets S.S. Basket, Drawer System, Oxygen Flow Meter. Flexible Gas Tubing, Nozzles Etc. Can Also Be Provided.



* Carrying Warranty at Rajkot Office only on site warranty not provided free of cost.

* Due to constant upgradation design, price and features are subject to change any time without prior notice.



WESTERN SURGICAL
EQUIPMENTS PRIVATE LIMITED

Western Surgical

96

* www.westernsurgical.in *
Email : westernsurgical@gmail.com
Ph. +91-281-2232871 / 2243159
Mo. +91 94280 03750/52/54

© Copyright Reserved by : Western Surgical - Rajkot.