

Model Popular (1+1) Std.



► Construction

Tubular rigid steel structure with detachable capacious drawer on smooth channels, powder coated for high corrosion resistance and aesthetic look.

► Monitor Tray

Top Monitor tray at eye level with a load capacity of 25 kg.

► Rotameters

In length 230 mm, rotating bobbin flowmeters calibrated in double / triple scale to ensure better accuracy and clarity in reading.

► Oxygen 100 cc/min to 8000 cc/min

► Nitrous Oxide 200 cc/min to 12000 cc/min.

► Vaporizer

1 No. Halothane Vaporizer & provision for other vaporizer. Emergency Press on block of oxygen Emergency oxygen cum outlet with pressure relief valve facility.

► Table Top

Table top made of superior quality stainless steel

► Cylinder Yoke

Pin index yoke with sliding steel clamping bars for easy handling each for O₂ and N₂O cylinders (5 Ltr water Capacity cyls.)

► Pressure Gauge

Superior quality pressure gauge of O₂ and N₂O panel mounted to read inside pressure of cyls.

► Regulators

Preset regulators for O₂ and N₂O with output pressure of 4.22 kg/cm². N₂O regulators are activated only when oxygen regulators are pressurised to minimize the risk of administering 100% N₂O.

► Metal Tubing

For easy mobility, the structure is provided on 4nos. large diameter castor wheels with ball bearing facility.

► Other facility

Provision to accommodate 2 ward type cylinders (Two 10 ltr. Water Capacity cylinders) of O₂ or N₂O cylinders.

► Provision for Circle Absorber fixing.

► Dimension

Height 1530 mm, Depth 540 mm, width 510 mm

► Weight

45 Kg (approx. without cyls.)

**Note: All models are available in Superior Quality S.S. at an extra cost.
Circle Absorber Both are Extra cost.**

* Carrying Warranty at Rajkot Office only on site warranty not provided free of cost.

* Due to constant upgradation design, price and features are subject to change any time without prior notice.



Western Surgical

WESTERN SURGICAL
EQUIPMENTS PRIVATE LIMITED

* www.westernsurgical.in *

Email : westernsurgical@gmail.com

Ph. +91-281-2232871 / 2243159

Mo. +91 94280 03750/52/54

© Copyright Reserved by : Western Surgical - Rajkot.