

ANAESTHESIA MACHINE (2+2)



- **Construction** : Tubular rigid steel structure with detachable capacious drawer on smooth channels, powder coated for high corrosion resistance and aesthetic look.
- **Monitor Tray** : Top Monitor tray at eye level with a load capacity of 25 kg.
- **Rotameters** : In length 230 mm, rotating bobbin flowmeters calibrated in double / triple scale to ensure better accuracy and clarity in reading.
 - Oxygen 100 cc/min to 8000 cc/min
 - Nitrous Oxide 200 cc/min to 12000 cc/min.
- **Vaporizer** : 1 No. Halothane Vaporizer & provision for other vaporizer.
- **Emergency Press on block of oxygen**
Emergency oxygen cum outlet with pressure relief valve facility.
- **Table Top** : Table top made of superior quality stainless steel

- **Cylinder Yoke** : Pin index yoke with sliding steel clamping bars for easy handling each for O2 and N2O cylinders (5 Ltr water Capacity cyls.)
- **Pressure Gauge** : Superior quality pressure gauge of O2 and N2O panel mounted to read inside pressure of cyls.
- **Regulators** : Preset regulators for O2 and N2O with output pressure of 4.22 kg/cm². N2O regulators are activated only when oxygen regulators are pressurised to minimize the risk of administering 100% N2O.
- **Metal Tubing** : For easy mobility. the structure is provided on 4nos. large diameter castor wheels with ball bearing facility.
- **Other facility** : Provision to accommodate 2 ward type cylinders (Two 10 ltr. Water Capacity cylinders) of O2 or N2O cylinders.
- **Provision for Circle Absorber fixing.**
- **Dimension** : Height 1530 mm, Depth 540 mm, width 510 mm
- **Weight** : 45 Kg (approx. without cyls.)

Note : All models are available in Superior Quality S.S. at an extra cost.

Circle Absorber Both are Extra cost.

* Carrying Warranty at Rajkot Office only on site warranty not provided free of cost.

* Due to constant upgradation design, price and features are subject to change any time without prior notice.



Western Surgical

WESTERN SURGICAL
EQUIPMENTS PRIVATE LIMITED

* www.westernsurgical.in *

Email : westernsurgical@gmail.com

Ph. +91-281-2232871 / 2243159

Mo. +91 94280 03750/52/54

© Copyright Reserved by : Western Surgical - Rajkot.