


Airway management

A range of airway management devices including a comprehensive choice of supraglottic airways, from the innovative i-gel to laryngeal mask options. The range also includes sterile and unsterile Guedel airways, endotracheal tubes, laryngoscope blades, filters, HMEs HMEFs and patient connections.

* **Airway devices**

* **Patient Connections**

* **Breathing filters, HMEs and HMEFs**



INTERSURGICAL i-gel[®] supraglottic airway

For use in anaesthesia and resuscitation



Code	8205	8204	8203	8225	8202	8215	8201
Box qty	25	25	25	10	10	10	10
Size	5	4	3	2.5	2	1.5	1
Weight (kg)	90+	50-90	30-60	25-35	10-25	5-12	2-5
Patient	large adult	medium adult	small adult	large paediatric	small paediatric	infant	neonate

The i-gel[®]

A truly unique, single use, supraglottic airway with a non-inflatable cuff, available in a full range of adult and paediatric sizes. i-gel incorporates a gastric channel to improve and enhance patient safety, an integral bite block to reduce the possibility of airway channel occlusion and buccal cavity stabiliser to aid insertion and avoid the potential for rotation. i-gel has been included in the new European Resuscitation Council Guidelines published in 2010.




www.i-gel.com

visit www.i-gel.com for videos, user guides, product information, FAQs and other downloads

Accurate and natural positioning

The i-gel accurately and naturally positions itself over the laryngeal framework, providing a reliable perilaryngeal seal without the need for an inflatable cuff.



* Carrying Warranty at Rajkot Office only on site warranty not provided free of cost.

* Due to constant upgradation design, price and features are subject to change any time without prior notice.



WESTERN SURGICAL
EQUIPMENTS PRIVATE LIMITED

Western Surgical

19

* www.westernsurgical.in *

Email : westernsurgical@gmail.com

Ph. +91-281-2232871 / 2243159

Mo. +91 94280 03750/52/54

© Copyright Reserved by : Western Surgical - Rajkot.