

A Fibre-Optic Laryngoscope Blade That Fits On Your Hook-on Handle !

FEATURES

- Hook-on Handle used with Conventional Blades can be used with Fibrelite Blades also. Hence no need to purchase any extra or special Handle.
- No Trauma causing Hot Bulb
- No possibility of loose bulb entering oral airway
- Cool high power light
- No Electrical Wires in the blade
- No Electrical Switches
- Full Stainless Steel construction
- No Flickering
- Over 3000 fibres bundled together ensures steady, bright & focused beam

The Handle shown in the picture is sold separately. There is no need to purchase any extra or special handle, as your existing hook-on handle can be used.



Macintosh Fibrelite Blades (English Profile)

Size	Use	Length	Code (Shiny)	Code (Matt)
1	Infant	92 mm	N/A	6509M
2	Child	100 mm	N/A	6507M
3	Adult	130 mm	N/A	6506M
4	Large Adult	155 mm	N/A	6509M

Miller Fibrelite Blades

Size	Use	Length	Code (Shiny)	Code (Matt)
0	Neonatal	25 mm	N/A	6516M
1	Infant	102 mm	N/A	6517M

Fibrelite Lamp
6571



Fibrelite Lamp Key
6584



Only Hook-on type Blades & Handles are marketed by us. Fibrelite Blades are only available in Matt finish. All Blades are made from Stainless Steel. All Handles are made from Chrome plated Brass & available in Shiny finish only.

5

Trupiti Fibrelite Blades with Flexible Articulating Tip

FEATURES

- Articulating tip provides precise control to elevate epiglottis
- Built in channel to help visualise and intubate with an unobstructed view
- Unique anti-seize spring action ensures long life & smooth operation
- Can be used on any standard hook on handle complying to ISO standards, hence no need to buy an additional handle separately
- An Innovative Indigenous Product that guarantees prompt after sales service locally with easily available genuine spares unlike in the case of imported counterparts
- Fully constructed out of Stainless Steel with matt finish

The Handle shown in the picture is sold separately. There is no need to purchase any extra or special handle, as your existing hook-on handle can be used.



Unbeatable
value for money
EASILY AFFORDABLE
even for the
"Just out of College"
Anaesthesiologist.



Trupiti Fibrelite Blades (available only in matt finish)

Size	Use	Length	Code
2	Child	100 mm	6522
3	Adult	130 mm	6523
4	Large Adult	155 mm	6524

Green System Fibre-Optic

Green System Trupiti Blades (Fibre-optic) with Flexible Articulating Tip

FEATURES

- Articulating tip provides precise control to elevate epiglottis
- Built in channel to help visualise and intubate with an unobstructed view
- Unique anti-seize spring action ensures long life & smooth operation
- An Innovative Product that guarantees prompt service with easily available genuine spares
- Fully constructed out of Stainless Steel with matt finish

The Handle shown in the picture is sold separately.



Trupiti Green System Blades Macintosh (American Profile)

Size	Use	Length	Code
2	Child	108 mm	6822
3	Adult	128 mm	6823
4	Large Adult	159 mm	6824

Trupiti Green System Blades Macintosh (English Profile)

Size	Use	Length	Code
2	Child	100 mm	6722
3	Adult	130 mm	6723
4	Large Adult	155 mm	6724

Standard
Handle
Accepts two size
"C" 1.5 V cells.



Stubby
Handle
Accepts two size
"AA" 1.5 V cells.



Penlite
Handle
Accepts two size
"AA" 1.5 V cells.



Howland
Lock



Green System Fibre-Optic

Green System

Macintosh Blades American Profile

The American profile with a slightly more abrupt curvature and a flange starting at the lock of the blade and continuing to approximately where the curvature changes direction then disappearing.

Size	Use	Length	Code
0	Neonatal	80 mm	6805
1	Infant	87 mm	6806
2	Child	108 mm	6808
2.5	Large Child	118 mm	6809
3	Adult	128 mm	6810
3.5	Medium Adult	144 mm	6811
4	Large Adult	159 mm	6812



Macintosh Blades English Profile

The Macintosh laryngoscope blade is one of the most popular blades in use today. The English profile has a smooth, gentle curve that extends from the base to the tip and a flange that starts at the lock and continues to the tip, narrowing gradually.

Size	Use	Length	Code
1	Infant	92 mm	6906
2	Child	100 mm	6907
3	Adult	130 mm	6908
4	Large Adult	155 mm	6909
5	Extra Large Adult	178 mm	6910



Miller Blades

One of the most popular of the straight blades. The tongue is straight with a slight upward curve near the tip. In cross section, the flange web and tongue form a C with the top flattened.

Size	Use	Length	Code
0	Neonatal	77 mm	6716
1	Infant	102 mm	6717
2	Child	153 mm	6718
3	Adult	195 mm	6719
4	Large Adult	205 mm	6720



* Carrying Warranty at Rajkot Office only on site warranty not provided free of cost.

* Due to constant upgradation design, price and features are subject to change any time without prior notice.



WESTERN SURGICAL
EQUIPMENTS PRIVATE LIMITED

Western Surgical

21

* www.westernsurgical.in *

Email : westernsurgical@gmail.com

Ph. +91-281-2232871 / 2243159

Mo. +91 94280 03750/52/54

© Copyright Reserved by : Western Surgical - Rajkot.