

Howland Lock

Most laryngoscope blades form a 90° angle with handle when ready for use. When a Howland lock is fitted between the handle and the blade by means of Hook-on connections on both ends, it forms a 45° angle. This change in the angle helps natural lifting action which simplifies natural exposure making intubation possible in the most difficult of cases.

The Handle & Blade shown in the picture are sold separately.



6175

Standard Handle
Accepts two size "C" 1.5 V cells.

6001



6003



Penlite Handle
Accepts two size "AA" 1.5 V cells.

6002



Stubby Handle

Useful during difficult intubation on a short neck, barrel chested or large bodied patient where the distance between the mouth & the chest restricts the use of a standard handle. Stubby accepts two size "AA" 1.5 volts cells.



15

Intubation Instruments

LMA Unique - Disposable

Manufactured by: The Laryngeal Mask Co. Ltd.

Used successfully on millions of patients worldwide, the LMA Unique can be used for a wide range of routine applications ranging from general anaesthesia to emergency use or as a resuscitation device.

- Sterile
- Muscle relaxants not required for mask insertion
- Smooth emergence from anaesthesia
- Minimal haemodynamic response
- Minimal soft tissue and dental trauma
- A secure, hands free airway
- Disposable
- Latex free



LMA UNIQUE is the trade mark of The Laryngeal Mask Company Limited.

Code	7714A	7714B	7714C	7714D	7714E	7714F	7714G
Size	1	1.5	2	2.5	3	4	5
Patient Selection Information	Neonate < 5kg	Infant 5-10kg	Child 10-20kg	Child 20-30kg	Child 30-50kg	Adult 50-70kg	Adult > 70kg
Maximum Inflation Volume	4ml	7ml	13ml	14ml	20ml	30ml	40ml

Laryngoscopes Conventional

Macintosh Blades (American Profile)

The American profile with a slightly more abrupt curvature and a flange starting at the lock of the blade and continuing to approximately where the curvature changes direction then disappearing.

Size	Use	Length	Code (Shiny)	Code (Matt)
0	Neonatal	80 mm	6705	6705M
0.5	Infant	87 mm	6706	6706M
1	Child	108 mm	6708	6708M
2	Large Child	118 mm	6709	6709M
3	Adult	128 mm	6710	6710M
3.5	Medium Adult	144 mm	6711	6711M
4	Large Adult	159 mm	6712	6712M



Macintosh Blades (English Profile)

The Macintosh laryngoscope blade is one of the most popular blades in use today. The English profile has a smooth, gentle curve that extends from the base to the tip and a flange that starts at the lock and continues to the tip, narrowing gradually.

Size	Use	Length	Code (Shiny)	Code (Matt)
00	Premature	70 mm	6004	6004M
0	Neonatal	80 mm	6005	6005M
1	Infant	92 mm	6006	6006M
2	Child	100 mm	6007	6007M
3	Adult	130 mm	6008	6008M
4	Large Adult	155 mm	6009	6009M
5	Extra Large Adult	178 mm	6010	6010M



Miller Blades

One of the most popular of the straight blades. The tongue is straight with a slight upward curve near the tip. In cross section, the flange web and tongue form a C with the top flattened.

Size	Use	Length	Code (Shiny)	Code (Matt)
00	Premature	57 mm	6015	6015M
0	Neonatal	75 mm	6016	6016M
0.5	Small Infant	88 mm	6014	6014M
1	Infant	102 mm	6017	6017M
2	Child	153 mm	6018	6018M
3	Adult	195 mm	6019	6019M
4	Large Adult	205 mm	6020	6020M



* Carrying Warranty at Rajkot Office only on site warranty not provided free of cost.

* Due to constant upgradation design, price and features are subject to change any time without prior notice.



WESTERN SURGICAL
EQUIPMENTS PRIVATE LIMITED

Western Surgical

23

* www.westernsurgical.in *

Email : westernsurgical@gmail.com

Ph. +91-281-2232871 / 2243159

Mo. +91 94280 03750/52/54

© Copyright Reserved by : Western Surgical - Rajkot.