

## Flexible Catheter Mounts



## 7042 - With Double Swivel All Plastic Disposable (APD)

Set comprises : 15cm X 11mm Corrugated Plastic Hose (2039), 15mm/22mm M/M Double Swivel (7005), 22mm Female Catheter Mount (9008)

## 7022 - With Double Swivel Plastic &amp; Rubber Reusable (PRR)

Set comprises : 15cm X 11mm Corrugated Rubber Hose (8102), 15mm/22mm M/M Double Swivel (7005), 22mm Female Catheter Mount (9008), 15mm Female Catheter Mount (9021F)

## 7044 - With Angle Mount All Plastic Disposable (APD)

Set comprises : 15cm X 11mm Corrugated Plastic Hose (2039), 15mm/22mm Angle Mount (9004), 22mm Female Catheter Mount (9008)

## 7024 - With Angle Mount Plastic &amp; Rubber Reusable (PRR)

Set comprises : 15cm X 11mm Corrugated Rubber Hose (8102), 15mm/22mm Angle Mount (9004), 22mm Female Catheter Mount (9008), 15mm Female Catheter Mount (9021F)



2707



2708



2701



2064



2060



2706



2062



2056



2063



## Snuggly Facemask

| Sizes         | All Silicone | All Black |
|---------------|--------------|-----------|
| 0 Infant      | 8300         | 8350      |
| 1 Small Child | 8301         | 8351      |
| 2 Child       | 8302         | 8352      |
| 3 Small Adult | 8303         | 8353      |
| 4 Adult       | 8304         | 8354      |
| 5 Large Adult | 8305         | 8355      |

A unique autoclavable snug facemask which moulds virtually into any shape giving unparalleled patient seal, has a soft contoured cuff for superior patient's comfort & easy to clean because of one piece construction. Silicone version allows clear view of patient lip colour / presence of vomitus / secretions.

## Clearhood

| Sizes         | All Silicone | All Black |
|---------------|--------------|-----------|
| 3 Small Adult | 8051         | 8044      |
| 4 Adult       | 8052         | 8045      |
| 5 Large Adult | 8053         | 8046      |

An ideal light weight transparent hooded facemask for resuscitation and anaesthesia. The strong clear polycarbonate hood allows for early recognition of vomitus, blood, foreign body and lets the anaesthetist monitor lip colour and intermittent misting. The supersoft cuff in silicone or black rubber gives a comfortable, least traumatic and safe seal. The cuff can be easily removed for cleaning/autoclaving.



## Circular Paedi

| Sizes          | All Silicone | All Black |
|----------------|--------------|-----------|
| 00 Premature   | 8014         | 8008      |
| 0 Infant       | 8015         | 8009      |
| 0A Small Child | 8016         | 8010      |
| 0B Child       | 8017         | 8011      |

One Piece design ensures gentle yet excellent fit. Available in silicone or black rubber versions.

## Rendell Baker Soucek

| Sizes         | All Silicone | All Black |
|---------------|--------------|-----------|
| 0 Premature   | 8036         | 8030      |
| 1 Infant      | 8037         | 8031      |
| 2 Small Child | 8038         | 8032      |
| 3 Child       | 8039         | 8033      |

The RBS mask cuts dead space to a quarter of anatomical facemasks with inflatable pad & hence give greater safety margin. Available in black & silicone versions. Silicone facemasks offer clear visualisation & patient's lip colour.

## Rubber Accessories



## Anatomical with Inflatable Pad

| Sizes         |  |      |
|---------------|--|------|
| 0 Infant      |  | 8056 |
| 1 Small Child |  | 8057 |
| 2 Child       |  | 8058 |
| 3 Small Adult |  | 8059 |
| 4 Adult       |  | 8060 |
| 5 Large Adult |  | 8061 |

Has a body made of rubber that can be widened or narrowed to fit the face with inflatable pad to ensure best fit & seal available.

## Head Harness - 8072

Several holes permit ventilation & holds both the Child & Adult heads. It consists of four rubber straps for attaching to the metal Mask Ring, which fits over the facemask.

## Mask Ring S.S. - 8071

Corrugated Breathing Tubes  
Silicone Rubber  
(Without Mounts)

| Sizes         |  |      |
|---------------|--|------|
| 10 cm x 11 mm |  | 8110 |
| 15 cm x 11 mm |  | 8111 |
| 30 cm x 11 mm |  | 8112 |
| 30 cm x 22 mm |  | 8113 |

Made from flexible light weight silicone rubber for extended life.

Corrugated Breathing Tubes  
Black Rubber  
(Without Mounts)

| Sizes          |  |      |
|----------------|--|------|
| 10 cm x 11 mm  |  | 8101 |
| 15 cm x 11 mm  |  | 8102 |
| 30 cm x 11 mm  |  | 8103 |
| 30 cm x 22 mm  |  | 8104 |
| 105 cm x 22 mm |  | 8107 |

Made from flexible light weight black rubber for extended life.

## Rubber Accessories

## Breathing Bags Black

| Sizes   | Neck       |       |
|---------|------------|-------|
| 0.5 ltr | 15 mm Bush | 8201B |
| 1.0 ltr | 22 mm Bush | 8202B |
| 1.5 ltr | 22 mm Bush | 8203B |
| 2.0 ltr | 22 mm Bush | 8204B |
| 3.0 ltr | 22 mm Bush | 8205B |

Designed for more accurate visual monitoring and increased tactile sensation. Dual purpose tail may be left intact for closed end use or cut to provide Bag Bleed Valve.

## Breathing Bags Green

| Sizes   | Neck       |       |
|---------|------------|-------|
| 0.5 ltr | 15 mm Bush | 8201G |
| 1.0 ltr | 22 mm Bush | 8202G |
| 1.5 ltr | 22 mm Bush | 8203G |
| 2.0 ltr | 22 mm Bush | 8204G |

Designed for more accurate visual monitoring and increased tactile sensation. Dual purpose tail may be left intact for closed end use or cut to provide Bag Bleed Valve.

## Bag Bleed Valve - 8221

Manufactured of light weight delrin with adjustable orifice. This valve inserts into the tail of a Breathing Bag.

## Rao's Tube Tree - 8120

Prevents drag and pull on mask or tracheal tube and helps prevent accidental extubation.

Note : Silicone Rubber is long lasting and resistant to ageing and can be autoclaved regularly. Natural rubber devices are autoclavable but this will reduce their life span.



\* Carrying Warranty at Rajkot Office only on site warranty not provided free of cost.

\* Due to constant upgradation design, price and features are subject to change any time without prior notice.



**WESTERN SURGICAL**  
EQUIPMENTS PRIVATE LIMITED

Western Surgical

26

\* [www.westernsurgical.in](http://www.westernsurgical.in) \*

Email : [westernsurgical@gmail.com](mailto:westernsurgical@gmail.com)

Ph. +91-281-2232871 / 2243159

Mo. +91 94280 03750/52/54

© Copyright Reserved by : Western Surgical - Rajkot.