



Jet Ventilator for Bronchoscopy

3026 - A ventilating system for rigid wall bronchoscopy and laryngeal microsurgery. Comprises : 2mtrs. High Pressure Tubing, Thumb Valve Control and three Luer Lock Needles.

Rao's Twin Gauge Oxygen Regulator

3027 - Specially introduced for use with the Jet Ventilator and other applications, it has an extra gauge to indicate the outlet pressure of the regulator and incorporates a safety valve. The outlet pressure can be adjusted with the help of the key to desired value depending upon use for Paediatric or Adult patients.



Spanners & Keys

3023 - Combination Key & Spanner - For fitting Regulators & F/A Valves onto Cylinders and for opening Bull Nose type Cylinders. Cannot be used to open Pin Index type Cylinders.

3028 - Multipurpose - For fitting flowmeters onto F/A Valves or Regulators, Hex Nuts onto Bull Nose type & Flush type Cylinders & for opening and closing both Bull Nose type as well as Pin Index type Cylinders.

3024 - Key for Pin Index Cylinders - For opening Pin Index type Cylinders & fitting Hex Nuts to the valves.

Twin Flowmeter Adaptor

3011 - For using two BPC flowmeters from one oxygen regulator.

PIPELINE ACCESSORIES



Pipeline Accessories



Terminal Units

3460 - Oxygen (White)
3462 - Air (Grey)
3461 - Nitrous Oxide (Blue)
3463 - Vacuum (Yellow)

Direct Probes

3450 - Oxygen (White)
3452 - Air (Grey)
3451 - Nitrous Oxide (Blue)
3453 - Vacuum (Yellow)

FEATURES

- Allows plugging of oxygen flowmeter & probes from front.
- Bears weight of the metal BPC flowmeter & humidifier.
- Push to insert and press-to-release mechanism for flowmeters & probes.
- Locating terminal unit pin & probe notch gives maximum stability to the flowmeter unit.
- Non-interchangeability between probes & terminal outlets for different services; possibility of error in making connection is avoided.
- Safety key plug used with Schrader type self sealing valves can also be used on the terminal units so that existing ancillary equipment can be inserted in both the terminal units as well as Schrader type sealing.



Single pipeline adaptor

3301 - Oxygen
Making & breaking the supply connection is in the same manner as described under "Self Sealing Valves"

Twin pipeline adaptor

3302 - Oxygen
Making & breaking the supply connection is in the same manner as described under "Self Sealing Valves"

Pipeline Accessories

Schrader Self Sealing Valves & Accessories

Illustration shows complete Oxygen Outlet Point (with Twin Pipeline Adaptor & BPC Flowmeters), as would be seen in actual use in wards.

Self sealing valves (Schrader type)

3303 - Oxygen
3305 - Air
3304 - Nitrous Oxide
3306 - Vacuum

The outlet points from which gases & vacuum are drawn are of self sealing type and have non-interchangeable outlet fittings. The supply is established simply by inserting the appropriate safety key plug into the orifice of the outlet point and pressing it home until it locks in position. To break the connection, the milled collar round the orifice is rotated a few degrees to the right, when the safety key plug is automatically ejected and the gas flow is cut off.

Safety key plugs (Schrader type)

3307 - Oxygen
3309 - Air
3308 - Nitrous Oxide
3310 - Vacuum

These are non-interchangeable like the self sealing valves and thus eliminate the possibility of errors in making connections.

Button cover

3316 - To keep the covers in position on the 'S' bracket.

'S' Bracket

3315 - For connecting the down drop pipes to self sealing valves.

Plastic covers for self sealing valves

3311 - Oxygen
3313 - Air
3312 - Nitrous Oxide
3314 - Vacuum
Coloured coded text for respective self sealing valves to safe guard against the possibility of administering the wrong gas.

* Carrying Warranty at Rajkot Office only on site warranty not provided free of cost.

* Due to constant upgradation design, price and features are subject to change any time without prior notice.



WESTERN SURGICAL
EQUIPMENTS PRIVATE LIMITED

Western Surgical

* www.westernsurgical.in *

Email : westernsurgical@gmail.com

Ph. +91-281-2232871 / 2243159

Mo. +91 94280 03750/52/54

© Copyright Reserved by : Western Surgical - Rajkot.