

ROTAMETER Hypoxia Guard



DIAMENSION - 3

Special Features :

- * Minimum 25% (+5%) Oxygen (1:3 mixed ratio between O₂ & N₂O) in surgery gas mixture.
- * 9" Rota meter Tube of O₂ & N₂O each also included in this unit.
- * Hypoxia Guard is a pneumatic device that ensure no. Nitrous Oxide Flow. if oxygen flow is lower than 0.5 liter/min approx.
- * Warranty : 12 Months from the date manufacturing

Circle Absorber Deluxe



CIRCLE Absorber Regular



Special Features :

SPECIFICATION

- * Double Chamber Adult Size
- * Easy cleaning - fully autoclavable.
- * Soda lime capacity 1800 gsm,
- * Transparent Light weight canister.
- * Completely leak proof.
- * Reversible Soda lime canister fitted to control head by means of a central tube which is screwed the head with a large plastic hand wheel.
- * Warranty : 12 month from the date manufacturing

Top Tech Circle Absorber



SINGLE NEW

FEATURES :

- * Passive gas scavenging facility
- * Adjusting Pressure Limit (APL) valve.
- * Interchangeable facility for single & double chamber (Canister) compatible for pediatric & adult surgery.
- * Soda Lime capacity for each chamber (Canister) 0.8 Kg.
- * Easy refiling Reversible soda lime Chamber (Canister).
- * International Standard 22 mm Connector with standard accessories.
- * Easy cleaning, fully autoclavable, light weight & transparent chamber (Canister)
- * All joints are leak Proof through taper mechanism.



DOUBLE NEW

* Carrying Warranty at Rajkot Office only on site warranty not provided free of cost.

* Due to constant upgradation design, price and features are subject to change any time without prior notice.



Western Surgical

WESTERN SURGICAL
EQUIPMENTS PRIVATE LIMITED

* www.westernsurgical.in *

Email : westernsurgical@gmail.com

Ph. +91-281-2232871 / 2243159

Mo. +91 94280 03750/52/54

© Copyright Reserved by : Western Surgical - Rajkot.