



Minic SM_001 M

Capacity: -550 mmHg $\pm$ 10 at 12-15 LPM
Noise Level: <50 dB A $\pm$ 3 Almost whispers
Pump type: Oil Free Piston Pump
Housing: Shock Proof FRP Cabinet
Jars : Polycarbonate 2x2.5 liters capacity with mechanical over flow system & ABS Lid
Filter: Reusable bacterial filter
Tubing: Non - Collapsible Suction Tubing
Power: 220/230V AC, 50Hz, 1Ph
Vacuum Gauge: 2.0 inch, 0-760 mmHg

M.R.P. 12500/-

Nano - AC/DC

SM_002 N AC/DC

Capacity: -580 mmHg $\pm$ 10 at 12-15 LPM
Backup : 120 mins
Pump type: Oil Free Piston Pump
Housing: Complete ABS (Plastic) Cabinet
Jars : Polycarbonate 2x2.0 liters capacity with mechanical over flow system & ABS Lid
Filter: Reusable bacterial filter
Tubing: Non - Collapsible Suction Tubing
Power: 220/230V AC, 50Hz, 1Ph & 12V DC
Vacuum Gauge: 2.0 inch, 0-760 mmHg

M.R.P. 21000/-



E-VAC

Capacity: -700mmHg $\pm$ 10 at 20-30 LPM
Noise Level: <50 dB A $\pm$ 3 Almost whispers
Pump type: Noiseless Rotary Vane Pump
Housing: Powder Coated MS Cabinet
Jars : Glass 2x2.0 liters capacity with mechanical over flow system & ABS Lid
Filter: Reusable bacterial filter
Tubing: Non - Collapsible Suction Tubing
Power: 220/230V AC, 50Hz, 1Ph
Vacuum Gauge: 2.5 inch, 0-760 mmHg

M.R.P. 21,500/-



* Carrying Warranty at Rajkot Office only on site warranty not provided free of cost.

* Due to constant upgradation design, price and features are subject to change any time without prior notice.



Western Surgical

WESTERN SURGICAL
EQUIPMENTS PRIVATE LIMITED

* www.westernsurgical.in *

Email : westernsurgical@gmail.com

Ph. +91-281-2232871 / 2243159

Mo. +91 94280 03750/52/54

© Copyright Reserved by : Western Surgical - Rajkot.