



Foot Suction Regular S.S. Top

SM_009 FSRT

M.R.P. 6000/-

Portable Twin glass of 1 Litr superior piston pump for creating vacuum instantly, base made from M.S. Sheet full powder coated, all aluminum parts are fully powder coated vacuum gauge 2½" gauge 2½" graduated in mmHg. (Also available in stainless steel base.)



Foot Suction Imported

SM_0011 FSI

Specification

Base : Power Coated Heavy MS Base
Capacity : -600mg Hg±10 at 200ml Per Stroke
Pump Type : Heavy Piston Pump
Jar : Polycarbonate 1500ml capacity
with over flow system & ABS Lid
Filter : Bacterial Filter Autoclavable/Reusable
Vacuum Gauge : 2.0 inch, 0-760mm Hg
Operated : Manual

M.R.P. 9500/-



Vacuum Extractor

SM_0012 VE

M.R.P. 15000/-

- Supreme Manual Vacuum Extractor
- Traction handle & Hand Suction pump
- Three Silicone cups (40,50,60mm)
- Glass jar fitted with vacuum gauge and control valve
- Optional:Stainless Steel instead of Silicone cups can be provided.
- Optional:Rubber tubes instead of silicone tubes can be provided.
- Stand for glass jar / • Silicone pressure tube



Silicone Cup SM_0013 SC

- 100% silicone medical product
- This vacuum suction cup could assist and shorten second stage of labor in occipital delivery
- Minimize the risk of Traumas, Cephalohematoma and Forceps delivery
- The soft cup can protect and decrease possible danger of damaging the fetus' head
- The handle and valve assemblies are manufactured of chrome plated brass and can be replaced

M.R.P. 3500/-

* Carrying Warranty at Rajkot Office only on site warranty not provided free of cost.

* Due to constant upgradation design, price and features are subject to change any time without prior notice.



Western Surgical

WESTERN SURGICAL
EQUIPMENTS PRIVATE LIMITED

* www.westernsurgical.in *

Email : westernsurgical@gmail.com

Ph. +91-281-2232871 / 2243159

Mo. +91 94280 03750/52/54

© Copyright Reserved by : Western Surgical - Rajkot.